

drive



aerocarbon

POWERCHAIR

Weighing just

12kg

Carbon Fibre
Technology

Warranty

5 Years

drivedevilbiss.co.uk

AeroCarbon POWERCHAIR

Carbon fibre technology delivers exceptional strength without the weight. As part of the Aero Carbon family, this sleek, modern powerchair is designed for effortless manoeuvrability, advanced comfort, and intuitive control with optional attendant control.

Specifications

Length	91cm (35.8")
Width	61cm (24")
Height	89cm (35")
Seat Depth	42cm (16.5")
Seat Width	44cm (17.3")
Seat Height	44cm (17.3")
Backrest Width	42cm (16.5")
Folded Size	73x31x76cm (LxWxH)
Ground Clearance	3.2cm (1.25")
Battery Supplied	24V 10Ah Lithium-ion
Motor Type	24V / 150W (DC Brushless)
Range*	9.5 miles
Maximum Speed*	4mph
Carry Weight (excl. battery)	12kg
Total Weight (incl. battery)	13.5kg (29.8lb)
Maximum User Weight	136kg (21.4 stone)

AeroCarbon Powerchair:
ULCFPCHAIR

Attendant Controller & Fixings:
12-EC8003Z-00

WARRANTY:
Frame - 5 years.
Electrical Components - 3 years.
Battery - 1 year.
Excludes wear & tear, & misuse.



drive

£

Contact us for more information

Our lightest
powerchair
weighs in at
just 12kg



Perfect for
transport, with
easy folding
and ideal for
car boots



Premium quality,
lightweight
carbon fibre
frame with a 5
year warranty



Air travel safe,
lightweight
lithium battery
with off board
charging



Optional
easy-to-fit
attendant
control
available

