

Wizemed

CE UK
CA

Carbon Fibre Power Wheelchair C100

Incredibly light,
just 16kg (2.5 st)



KEY FEATURES

GRADIENT COLOURS

Your journey, your shade.

TRIM COLOURS



Ocean
Blue



Bordeaux
Red



Forest
Green



Sunset Orange



Aurora Green



Dusky Lavender



Carbon Black

KEY FEATURES

LIGHTWEIGHT
AND MANOEUVRABLE

Full carbon fibre frame

Lightweight, high-strength, and corrosion-resistant.



Elastic backrest

Elastic breathable backrest, Ergonomic design for more comfortable sitting.



Advanced Control

Advanced controller with an electromagnetic brake for immediate stopping.

Electromagnetic Brake

Stops as soon as you let go, whether uphill or downhill.



Passed SGS Tests

Personalised Customisation

Multiple colour options to make your electric wheelchair stand out.



Aluminium Alloy Wheel Hub



Patented Three-Position Anti-Tip Wheels

TECHNICAL DATA

SPECIFICATIONS

Body Weight without battery	16 kg / 2.52 st / 35 lbs
Overall Dimensions(L*W*H)	950 × 580 × 900 mm / 37 × 23 × 35"
Folded Dimensions(L*W*H)	580 × 360 × 730 mm / 23 × 14 × 29"
Seat Width	440 mm / 17.3"
Seat Height(Cushion Removed)	460 mm / 18.1"
Seat Depth	420 mm / 16.5"
Maximum Speed	6 km/h / 3.7 mph

Range	15 km / 9.3 mi
Climbing Ability	0-8°
Turning Radius	800 mm / 31.5"
Maximum Load	136 kg / 21.4 st / 300 lbs
Motor Power	250w × 2 Brushless
Batteries	24 V / 10.4 Ah, Lithium-ion (Expandable)
Frame Material	Carbon Fiber

Copyright © 2024 Wizemed. All rights reserved.
Wizemed reserves the right to make changes to the design, specifications, and models without prior notice.

Wizemed B.V.
Suite 4.21-02 (4th floor)
Evert van den Beekstraat 354
1118 CZ Schiphol
The Netherlands

Wizemed Ltd
Regus Birmingham NEC
4200 Waterside Centre
Birmingham Business Park
Solihiull, B37 7YN

TEL: +447398458247
E-Mail: howie.wong@wizemed.com
www.wizemed.co.uk



Scan the QR
code to visit
Wizemed